

# Provo.X

**Provocation system**

Safe and hygienic challenge testing



# OVERVIEW

GANSHORN's Provo.X is an aerosol dosimeter that is able to do specific and unspecific provocations within seconds – no matter where the measurements are carried out. Perform inhalational provocation tests quickly and with reproducible results. The PC based system is designed for accurate and safe direct

challenge tests. Various protocols and sequences can be easily adapted. All accessories with patient contact are designed as disposable items, which 100 % excludes possible cross-contamination. Eliminating the disinfection steps save time and money.



**Configurable: precise and adjustable**



**Hygienic disposable solution**



**Provides standardized protocols and offers individual series**



**Flexible choice of control module**



**Powerful and quiet compressor**




**Stores up to six provocation series**

The measurement of hypersensitivity or hyper-responsiveness of the bronchial airway is necessary if asthma is suspected. With the help of bronchial provocation, certain lung diseases can be excluded. Specifically, a bronchial provocation test is indicated in the following cases:

- If there is a need to rule out hyper-reactivity as a diagnosis (e.g. asthma)
- Contradiction between asthmatic complaints and normal basic pulmonary function
- When results of provocation tests indicate expected therapeutic or preventive consequences

# LFX - SOFTWARE PLATFORM


The LFX software is GANSHORN's user-friendly interface. For provocation tests, LFX provides selectable standardized provocation steps for the management of controlled measurements. In addition, user defined provocation rows can be individually configured as required.



Select sequence name to display substance, residence, dose, etc., for every stage of the sequence

Select control module for type of measurement

The desired medication can be added, the dose is adjusted precisely by adapting the breaths and application time.




Adjust the residence time here, if needed. It will be displayed as countdown after the step.

Change the treatment time and/or the number of breaths in order to alter the delivered dose.

Add a new step, or adjust an existing one in a few seconds

Completed step  
Indicates next step to be taken  
Provocation substance, dose and treatment time and tabular results for each step



PD-20 = Provocation dose at which FEV1 for example, decreases by 20% (default) from the pre or base 100% reference. When FEV1 falls below PD-20 line, provocation should be stopped

Graphical comparison of bronchial challenges for each step


# Measurement principle

An FVC test is typically used for reference and post-provocation assessment. Other possible measurements for a provocation sequence include:

- Offline parameter (manual input)
- No measurement (e.g., biomarkers)
- Body plethysmography

Each control module allows for specific measurements and can be used to determine the PD reference with percentage differences from the base or pre-measurements.


## Reference measurement



**A**

## Inhalation of medication

The patient inhales an aerosolized provocation substance using a nebulizer, taking slow, deep breaths. The substance is nebulized during the inspiratory cycle for a set bolus time. If inhalation is shorter than bolus time, the remaining time is added to the next breath. LFX monitors breathing and administers the aerosol only during inhalation, optimizing delivery and minimizing waste. The total dose per stage is automatically adjusted for short breathing cycles. As nebulization progresses, water vapor increases the concentration of the test substance, so the volume of the nebulizer solution should be increased accordingly.




**B**

## Control test

After a set residence time (1-20 min), indicated by a countdown, the user undergoes a body plethysmography, spirometry or oscillometry test. Then, the next provocation stage can begin. This cycle of inhalation and control tests repeats until all stages are completed or the patient opts to stop.

For each control module, measurements can be defined and used for the PD reference with percentage differences from the base, such as a 20% decrease in FEV1 or a 100% increase in sRaw.



**C**

## Interpretation of results

If after the last inhalation the FEV1 or alternative control parameter doesn't fall below the threshold, the provocation test is regarded negative and asthma can be ruled out. If the threshold (PD20 in

spirometry) is reached, the provocation should be stopped and the PD20 value of this positive provocation test is reported.

# SET-UPS

## With SpiroScout

Combination with spirometer



## In PowerCube Diffusion+

For combined use with diffusion system



## Stand alone device

System with arm and clamp



## PowerCube Body+

For integrated use in bodybox






# FEATURES




## Threshold-controlled aerosol delivery

Provo.X administers the aerosol dose in a breath-synchronous manner during the inspiration phase, ensuring that the aerosol reaches deeply into the airways. That improves the reproducibility of results.



## User-Friendly Software

GANSHORN's software LFX includes pre-installed protocols, such as ATS and DGP, and allows for the creation and storage of custom protocols with up to 14 stages.



## Adjustable Dosage

The dosage can be adjusted by changing the concentration of the provocation liquid or the atomization time, providing precise control over the amount of aerosol delivered at each stage.

Step	Substance	Concentration	Breaths	Dose	Treatment time [s]
Pre	B				
1	NaCl	0.9 %	4	129 µg	0.49
2	Methacholin	3.3 mg/ml	3	15.8 µg	0.219
3	Methacholin	3.3 mg/ml	6	31.7 µg	0.219
4	Methacholin	3.3 mg/ml	10	63.4 µg	0.263
5	Methacholin	3.3 mg/ml	10	128 µg	0.53
6	Methacholin	3.3 mg/ml	14	240 µg	0.71
Lysis	Berotec		2	strokes	



## Hygienic Design

Components which can be contaminated are single use. They must not be sterilized what saves time and makes Provo.X one of a few devices with a perfect hygienic solution. This ensures maximum safety for patients.



# CONSTRUCTION

Due to the perfect user guidance and simple handling of the Provo.X, you can ease the daily work of your staff. The well-structured display of the Provo.X software LFX continuously provides precise information on the status of

the examination. The backside-ventilated special valve located directly at the atomizer and the large compressed air tank guarantee an optimal droplet size, which remains constant over the entire bolus time.

Mouth and T-piece



Provocation substance reservoir (single use)

Provo.X unit - communicates with the computer (LFX) via Bluetooth

Breath detection via pressure sensor

Compressed air in (e.g. from compressor)

The powerful compressor with a three-liter compressed air tank guarantees a stable atomizing pressure and a low noise level (operating time of only approx. 1-3 minutes per day).



# WHY GANSHORN?

For 40 years GANSHORN has been manufacturing a complete state-of-the-art portfolio of pulmonary function testing systems for spirometry, body plethysmography, diffusion, bronchial provocation and cardiopulmonary stress testing. With its technological innovations, the company has been a leader in the diagnostics market since 1982. Many of these

are now perceived as gold standards. In order to meet our high quality standards, it is important to us that all key components are made in Germany. Our devices are created in modern processes in Bavaria, from the initial idea to distribution. In the meantime GANSHORN is represented worldwide, with strong markets in Europe, Asia, North and South America.



## PowerCube Body+

Body Plethysmography



## SpiroScout

Spirometry



## PowerCube Diffusion+

Diffusion measurement



## Provo.X

Provocation testing



## PowerCube Ergo

Cardiopulmonary exercise testing (CPET)



## Vivatmo pro

FeNO monitoring



## tremoflo®

Airwave oscillometry



## EucapSys

EVH provocation



## AltiTrainer

Hypoxic challenge testing, hypoxia training

**GANSHORN Medizin Electronic GmbH**  
Industriestr. 6-8 | 97618 Niederlauer, Germany

✉ sales@ganshorn.de

☎ +49 9771 6222 0

🌐 www.ganshorn.de



019950157 – Rev 3.0



The model shown may also include optional equipment which is not within the standard scope of supply. Design, equipment, and contents are subject to change without notice, as are misprints and errors.